## Parcus Nitinol Guide Wire

### Important Product Information

Produktinformationsblatt

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## Parcus Nitinol Guide Wire

### 1. Indications:

The Parcus Nitinol Guide Wire is used to aid in the placement of cannulated interference screws.

### 2. Warnings:

- A. This product is intended for use by or on the order of a physician.
- B. This product is designed for single use only. Devices that are reused can fatigue and break.
- C. This product is not intended for use as an implant.

### 3. Material:

This device is manufactured from Nitinol per ASTM 2063. The materials used in the manufacture of this device that are intended to be placed inside the body are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

### 4. Sterilization:

This device is not provided sterile. The following table provides recommended minimum sterilization parameters that have been validated by Parcus Medical to provide a  $10^{-6}$  sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time
Gravity	132°C / 270°F	15 minutes	30 minutes <sup>1</sup>
Pre-vacuum	132°C / 270°F	4 minutes	

<sup>1</sup> Drying times vary according to load size and should be increased for larger loads.

### 5. Directions for use:

- A. Position graft in tunnel or socket in the preferred manner.
- B. Place guide wire in desired position.
- C. Place cannulated interference screw over guide wire.
- D. Engage cannulated screw driver with screw and turn screw into desired position.
- E. Remove driver.
- F. Remove and dispose of wire.

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