## Parcus Tibial Plug

### Important Product Information

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## Parcus Tibial Plug

### 1. Indications:

Parcus Tibial Plugs are used during ACL/PCL ligament repair or reconstruction procedures to limit the loss of irrigation fluid from the knee joint. The Parcus Tibial Plugs are supplied in two configurations; solid and with a 2.4mm cannulation. Both designs are tapered to accommodate tunnels from 7mm to 12mm in diameter.

### 2. Warnings:

- A. This product is intended for use by or on the order of a physician.
- B. Do not use this product if it shows any visible signs of damage.
- C. This device is not intended to be used to dilate the tibial tunnel. Care should be exercised when placing the plug to prevent fracture of the proximal tibia.
- D. This device is not intended to be used as an implant.

### 3. Material:

This device is comprised of PEEK, is radiolucent, and may be visualized through the use of ultrasound or MRI.

### 4. Inspection:

- A. Inspect the device for damage at all stages of handling.
- B. If damage is detected, consult the manufacturer for guidance.

### 5. Cleaning:

- A. Immediate rinsing and cleaning after use with an enzymatic detergent will effectively remove and prevent drying of adherent blood, tissue, etc.
- B. Scrub device with a soft brush, paying special attention to areas where debris might accumulate. Always avoid any harsh materials that can scratch or mar the surface of the device.
- C. Rinse the device thoroughly with water following the cleaning process.

### 6. Sterilization:

This product is a non-sterile item that must be cleaned and sterilized prior to use. The following table provides recommended minimum sterilization parameters that have been validated by Parcus Medical to provide a  $10^{-6}$  sterility assurance level (SAL):

| Cycle Type | Minimum Temperature | Minimum Exposure Time (Wrapped) | Minimum Dry Time        |
|------------|---------------------|---------------------------------|-------------------------|
| Gravity    | 132°C / 270°F       | 15 minutes                      | 30 minutes <sup>1</sup> |
| Pre-vacuum | 132°C / 270°F       | 4 minutes                       |                         |

<sup>1</sup> Drying times vary according to load size and should be increased for larger loads

### 7. Directions for use:

- A. Create a tibial tunnel that is 7mm to 12mm in diameter.
- B. If a guide wire is present in the tunnel, select the cannulated version of the Parcus Tibial Plug. The cannulation is 2.4mm in diameter. If no guide wire is present, the non-cannulated Plug can be used.
- C. Gently press the tapered end of the Plug into the tunnel until an interference fit, sufficient to keep the Plug in place, is achieved.
- D. Remove the Plug to pass instruments, graft material and suture.
- E. Replace the Plug, described in Step C, as necessary.

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