

*Surgical Innovation | Customer Driven*

**SLiK Fix** is a screw ideally suited for tenodesis procedures. The design of **SLiK Fix** screws provides procedural advantages compared to other fixation methods.

The **Tendon Washer and Trocar Tip** allow one to engage and advance the tendon to the appropriate depth into the bone socket and assess tension of the tendon prior to inserting the screw. The **Tendon Washer** provides an increased surface area to prevent tendon pull-back. The **SLiK Fix** driver features an **Activation Button** that enables a simple, one-handed technique to advance the fixation screw.

Recommended for use in both large and small-joint repairs.

### Key Features

- Anchor diameter options: 7.0mm, 8.0mm, 9.0mm
- Tendon Washer diameter: 5.0mm
- Tendon Washer/Trocar Tip on inserter to better manipulate the tendon
- Activation Button on driver for easy screw insertion
- Material solutions:
  - Carbon fiber-reinforced PEEK-OPTIMA from Invibio
  - PEEK-OPTIMA Natural from Invibio

**SLiK... because it is!**



### SLiK Fix PEEK CF Screws

Part #	Description
10995	7.0mm x 15mm <b>SLiK Fix</b> PEEK CF Screw
10996	8.0mm x 15mm <b>SLiK Fix</b> PEEK CF Screw
10997	9.0mm x 15mm <b>SLiK Fix</b> PEEK CF Screw

### SLiK Fix PEEK Screws

Part #	Description
10952	7.0mm x 15mm <b>SLiK Fix</b> PEEK Screw
10953	8.0mm x 15mm <b>SLiK Fix</b> PEEK Screw
10954	9.0mm x 15mm <b>SLiK Fix</b> PEEK Screw

### Required Instrumentation for SLiK Fix Screws

Part #	Description
11016	7mm <b>SLiK Fix</b> Headed Reamer, Sterile, Single-Use
11017	8mm <b>SLiK Fix</b> Headed Reamer, Sterile, Single-Use
11018	9mm <b>SLiK Fix</b> Headed Reamer, Sterile, Single-Use
11019	10mm <b>SLiK Fix</b> Headed Reamer, Sterile, Single-Use